ESACROM S.R.L.
VIA SERRAGLIO19 40026 IMOLA (BOLOGNA)
TEL 0542/643527 FAX 0542/647854
Web: www.esacrom.com E-Mail: esacrom@iol.it
www.esacrom.it

Esacrom srl

Manuale operativo Pulitrice

Versione 1/1 del 05/05/00

Brush 400

ESACROM S.R.L.
VIA SERRAGLIO19 40026 IMOLA (BOLOGNA)
TEL 0542/643527 FAX 0542/647854
Web: www.esacrom.com E-Mail: esacrom@iol.it
www.esacrom.it

Indice

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CAUTION

IMPORTANT!!!

ESACROM has compiled this booklet to help you handle and operate the equipment better.

ESACROM cannot accept any responsibility for damage resulting from incorrect use or modifications made to the equipment.

A copy of this booklet accompanies each machine. To ensure correct machine operation, please make sure you follow carefully the instructions provided.

The machine must never be used for purposes other than those indicated in the operating instructions and leaflet.

Non-compliance with the instructions contained in this booklet will result in the immediate termination of the guarantee.

On all requests, always indicate the date of purchase, the machine model and the serial number shown on the plate located on the back of the machine. Removal of this plate will result in the immediate termination of the guarantee.

All the technical details have to be intended as pure indicative information. The company Silfradent declares the right to modify the present manual without any compulsory previous advice to the customer.

INTRODUCTION

This manual refers to Brush 400 Polishing machine with double speed and a bench control.

ESACROM congratulates and thanks you for your choice.

The following instructions will help you get to know and use the machine better.

Always make sure the certificate of guarantee is attached to the machine and that this has been completed and stamped by the stockists. A copy of this guarantee should remain with the stockists and with ESACROM.

Before unpacking the machine, use the following check list to ensure no component parts are missing.

After carefully removing all parts from their packing, place these on a table to ensure they have not been damaged during transport. If any parts are found to be damaged, contact your stockists immediately. Also make sure that all packing residue is removed from the machine.

Now follow the installation instructions in this booklet.

EQUIPMENT

- CE Conformity
- Use and maintenance manual
- Warranty Certification
- Spindle for the right side
- Spindle for the left side
- Key to dismantle the spindles
- On request the right spindle to hold the disc

ESACROM S.R.L
VIA SERRAGLIO19 40026 IMOLA (BOLOGNA)
TEL 0542/643527 FAX 0542/647854
Web:www.esacrom.com E-Mail:esacrom@iol.it
www.esacrom.it

PRECAUTION

IMPORTANT - The Polishing machine is a dental laboratory instrument.

IMPORTANT - No other possible uses exist for this instrument outside the dental sphere.

The instrument must be used by qualified personnel who are familiar with the uses for which the instrument is intended.

The end user must indeed :

- Install the instrument on a plane and stable surface.
- Stop on working when you realize a damage of the feed cable or of the plug .
- wear always glasses or protection mask.

This instrument creates a circular movement at a high velocity, as consequence it may also be dangerous for the user and the persons behind him.

Don't never try to block this rotation with foreign bodies or with the hands.

Indeed you must consider that this is an equipment feed with 220 tension and you are recommended to take all the cautions of the case.

DESCRIPTION

The Polishing machine is used for the cleaning, templating and burnishing plaster models. The velocity is created by a protected motor which doesn't require any particular maintenance.

The wide base enables a high equilibrium ; it also has a self cooling system and four small rubber supports.

The motor can reach a rotation variable from 1400 rpm to 2800 rpm.

In case there is a sudden lack of current, the machine doesn't start automatically when the current is functioning again, it is in fact indispensable to repeat the switching on operation (CE Normative).

By pushing the general switch the two spindles start working at the same time.

The special plan of the instrument is intended to prevent any type of interference between the work zone and the motor.

INSTALLATION

Though during packing and preparation for shipment, all precautions have been taken to prevent damage during transport, it is important to check the condition of the machine after removing it from the packing. In the event of any damage, immediately notify the supplier.

It is best not to throw away the packaging as this may be useful if the machine has to be sent back for checking or overhauling.

IMPORTANT - Connect the feed lead to a power socket. Always ensure the power ratings are the same as those shown on the plate at the back of the machine.

Make sure the power socket is correctly earthed.

ESACROM S.R.L
VIA SERRAGLIO19 40026 IMOLA (BOLOGNA)
TEL 0542/643527 FAX 0542/647854
Web:www.esacrom.com E-Mail:esacrom@iol.it
www.esacrom.it

MACHINE SET-UP

The machine is sold ready to be used : it is sufficient to connect the feed cable to the machine and to set the spindles.

It is advisable to use the polishing machine with splash guards (Art. 906).

USE AND MAINTENANCE

The polishing machine can be used for all kinds of dental laboratory operations related to the polishing procedures.

The user can choose between the two shown velocities through the switch on the right side.

The minimum velocity is recommended for resins, while the maximum one for metals.

MAINTENANCE PROGRAMME

Concerning the maintenance of this machine, it is necessary to make a preliminary remark : inside there are bearings already lubricated, **as consequence there is no need of any additional lubricant.**

Indeed :

- Check there are no cuts, incisions etc, on the feed cable.
- Check if the brushes are efficient.

MAIN TECHNICAL PROBLEMS

A. If the Vibrator is out of order, what do you have to check ?

- 1) Check if the plug is well connected.
- 2) Check if you executed in the correct way the starting operation.
- 3) Probably the switch is faulty. We advise you to get in touch with a specialized technician.

B) What to do if during the use you hear a strange noise ?

- 1) Check if the noise is provoked by the right or the left bearing. If you realize this problem, we kindly ask you to get in touch with our assistance department.

SAFETY INSTRUCTIONS

- 1) Don't make any modification on the electrical system, in case of problems get in touch with a technician.
- 2) While you are cleaning the machine, remember to disconnect the current.
- 3) After the replacement of the brushes check carefully the instrument performance.
- 4) Follow the CE Normative for the end user security.

CAUTION

If there is the need to repair the instrument , we recommend you to apply to specialized centers. In case of replacements we recommend to use original spare parts (please make reference to the herewith enclosed exploded view).

ESACROM S.R.L
VIA SERRAGLIO19 40026 IMOLA (BOLOGNA)
TEL 0542/643527 FAX 0542/647854
Web:www.esacrom.com E-Mail:esacrom@iol.it
www.esacrom.it

TECHNICAL FEATURES

Voltage	230 V. +/- 15% (110 V on request)
Frequency	50/60 Hz
Absorbed power	350 Watt
Minimum velocity	1400 rpm
Maximum velocity	2800 rpm
Dimensions (L x W x H)	46x19x25 cm
Spindles height	17 cm.
Weight	11.35 Kg

SPARES

Right spindle	code	
Left spindle	code	
Condenser	code	
Switch to change the velocity	code	
General switch	code	
Bearings (left & right)	code	